

# Algorithm for the Resolution of Discrepancies of the ABO System

Presented by **Dra. Denise Harmening**  
**Director of Instruction & Customer Learning**  
**American Red Cross, Heritage Division**


**Professor, Department of Medical & Research Technology**  
**University of Maryland School of Medicine**  
**Baltimore, Maryland USA**

---

The need is constant.  
The gratification is instant.  
Give blood.™



# Who am I?



## Perspectives

## Historical

- 1901-Karl Landsteiner drew blood from himself and five co-workers
- separated the cells and serum; mixed each cell sample with each serum
- first to perform forward and reverse grouping

The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening



## ABO GROUPING: RECIPROCAL RELATION

- **FORWARD GROUPING**
- **REVERSE GROUPING**
- **EXPECTED REACTIVITY: 3+ or 4+**

# ABO ANTIBODIES

- **NATURALLY OCCURRING (expected)**
- **HIGH TITERED**
- **TYPICALLY IgM (small quantities of IgG may be present)**
  - **Exception: Group O which has a unique IgG Anti-A,B**
- **CANNOT CROSS PLACENTA (EXCEPT ANTI-A,B and IgG ANTI-A OR ANTI-B)**
- **BINDS & ACTIVATES COMPLEMENT**
- **CAPABLE OF CAUSING INTRAVASCULAR HEMOLYSIS**

# GROUP O INDIVIDUALS

- **Produce anti-A,B**
  - Typically IgG
  - Separate entity
  - Importance
    - Intravascular and extravascular transfusion reaction-HTR
    - May cross placenta--ABO HDN

# ABO GROUPING

## FORWARD

## REVERSE

<b>Anti-A Reagent</b>	<b>Anti-B Reagent</b>	<b>A1 cells Reagent</b>	<b>B cells Reagent</b>	<b>Interpret</b>
<b>4+</b>	<b>0</b>	<b>0</b>	<b>4+</b>	<b>A</b>
<b>0</b>	<b>4+</b>	<b>4+</b>	<b>0</b>	<b>B</b>
<b>4+</b>	<b>4+</b>	<b>0</b>	<b>0</b>	<b>AB</b>
<b>0</b>	<b>0</b>	<b>4+</b>	<b>4+</b>	<b>0</b>

The need is constant.  
The gratification is instant.  
Give blood.™


© D. Harmening



# Exception: Cord Blood / Neonates

- **Cord blood**
  - Wash cells before testing
  - Serum testing not commonly performed
- **Neonatal sample**
  - Serum testing not commonly performed

# Tube Testing: The Gold Standard




The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening



3+

SEVERAL LARGE  
AGGLUTINATES-CLEAR  
BACKGROUND




SCORE: 9

4+


ONE SOLID  
AGGLUTINATE



SCORE: 11



NEGATIVE: NO AGGREGATES



The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening




## Gel Technology


**How are reactions graded in  
blood banking?**

# GEL TECHNOLOGY


## ID-MTS REACTION GRADING CHART




A **negative** reaction is characterized by unagglutinated red cells forming a well-delineated pellet at the bottom of the microtube.




A **1+** reaction is characterized by red cell agglutinates predominantly observed in the lower half of the gel column. Unagglutinated cells form a pellet in the bottom of the microtube.



A **2+** reaction is characterized by red cell agglutinates dispersed throughout the length of the gel column. Few agglutinates may be observed in the bottom of the microtube.



A **3+** reaction is characterized by the majority of red cell agglutinates trapped in the upper half of the gel column.



A **4+** reaction is characterized by a solid band of red cell agglutinates on top of the gel. A few agglutinates may filter into the gel but remain near the predominant band.



A **mixed cell** reaction is characterized by a band of red cell agglutinates on top of the gel, accompanied by a pellet of unagglutinated cells at the bottom of the microtube.

The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening



# ABO SYSTEM: A QUICK REVIEW

## INHERITANCE CODOMINANT

- FOLLOWS MENDEL'S LAW OF INHERITANCE
- INDIVIDUALS INHERIT ONE GENE FROM EACH PARENT
- TWO GENES DETERMINE TYPE
- O GENE IS AN AMORPH
- PHENOTYPE IS DESCRIBED AS A, B, AB, OR O BLOOD GROUPS

# AB x AB MATING

GENES	<i>A</i>	<i>B</i>
<i>A</i>	<i>AA</i>	<i>AB</i>
<i>B</i>	<i>AB</i>	<i>BB</i>

The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening



# EXCEPTION: RARE *Cis-AB*

- Inheritance of both *AB* genes from one parent carried on one chromosome and an *O* gene inherited from the other parent
- Offspring inherits three *ABO* genes
- B antigen reacts weaker with anti-B
- Serum contains anti-B which will react with normal B cells, not *cis-AB*

# ABH ANTIGENS

- Referred to as glycolipids on Rbc membrane
- Referred to as glycoproteins in secretions
- Inherited genes produce specific **GLYCOLSYLTRANSFERASES** that add sugars to basic precursor substances
- Develop as early as 37 days of gestation, but not fully developed at birth (Newborn antigens react weaker until fully developed at 6-18 months of age)
- Present on lymphocytes, platelets, kidney, epithelium tissue, etc.

# RBC ANTIGEN EXPRESSION

- ABO Genes--chromosome 9
  - A and B encode, O is amorphic
- Hh Genes--chromosome 19
  - H encodes, h is amorphic

**GLYCOSYLTRANSFERASES AND IMMUNODOMINANT SUGARS RESPONSIBLE FOR H, A, AND B ANTIGEN SPECIFICITIES**


<b>Gene</b>	<b>Glycosyltransferase</b>	<b>Immunodominant Sugar</b>	<b>Antigen</b>
<b><i>H</i></b>	<b><math>\alpha</math>-2-L-fucosyltransferase</b>	<b>L-fucose</b>	<b>H</b>
<b><i>A</i></b>	<b><math>\alpha</math>-3-N-acetylgalactosaminyltransferase</b>	<b>N-acetyl-D-galactosamine</b>	<b>A</b>
<b><i>B</i></b>	<b><math>\alpha</math>-3-D-galactosyltransferase</b>	<b>D-galactose</b>	<b>B</b>

The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening




# PARAGLOBOSIDE



The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening

# H ANTIGEN



The need is constant.  
The gratification is instant.  
Give blood.™


© D. Harmening

# H Antigen

- 99.99% frequency in population
- *hh* - Bombay phenotype
- ABO antigens cannot be expressed if the H gene was not inherited
- ABO expression is dependent upon H inheritance

# B ANTIGEN


D-Galactose: Immunodominant Sugar responsible for "B" specificity



The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening

# A ANTIGEN



The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening

# RBC H ANTIGEN EXPRESSION

**O > A<sub>2</sub> > B > A<sub>2</sub>B > A<sub>1</sub> > A<sub>1</sub>B**

**most H -----> least H**

# What are the frequencies of ABO in the US Population?



The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening

# THE ABO BLOOD GROUP SYSTEM

**Table 5-3. ABO Phenotype Frequencies in U.S. Populations**

Phenotype	White (%)	Black (%)	Mexican (%)	Asian (%)
O	45	49	56	43
A <sub>1</sub>	33	19	22	27
A <sub>2</sub>	8	8	6	Rare
B	10	19	13	25
A <sub>1</sub> B	3	3	4	5
A <sub>2</sub> B	1	1	Rare	Rare

Garratty G, Glynn SA, McEntire R : ABO and Rh (D) phenotype frequencies of different racial/ethnic groups in the United States. Transfusion 44:703-706, 2004.

The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening



# DISCREPANCIES

- **No reciprocal relationship between the forward and reverse grouping**
- **You may have missing, extra, or weak Rx due to:**
  - **Technical errors**
  - **Unexpected reactions (reverse group- problems with serum)**
  - **Unexpected reactions (forward group- problems with red cells)**
  - **Unexpected reactions (both forward and reverse)**

# Common Sources of Technical Errors Resulting in ABO Discrepancies

- **Inadequate identification of blood specimens, test tubes, or slides**
- **Cells suspension either too heavy or too light**
- **Clerical errors**
- **A mix-up in samples**
- **Missed observation of hemolysis**

# Common Sources of Technical Errors Resulting in ABO Discrepancies

- **Failure to add reagents**
- **Failure to follow manufacturer's instructions**
- **Uncalibrated centrifuge**
- **Under centrifugation**
- **Over centrifugation**
- **Contaminated reagents**
- **Warming during centrifugation**

# ALGORITHM/ GENERAL GUIDELINES FOR RESOLUTION OF ABO DISCREPANCIES

- ALWAYS WASH PATIENT RBC SUSPENSION WITH NORMAL SALINE AND
- ALWAYS REPEAT THE TEST
- ALWAYS OBTAIN THE PATIENT'S HISTORY, AGE, DIAGNOSIS, TRANSFUSION HISTORY AND MEDICATIONS, IF POSSIBLE

## IF PROBLEM PERSISTS:

- TEST PATIENT'S CELLS WITH AVAILABLE LECTINS AND ANTI-A,B REAGENT
- TEST PATIENT'S SERUM WITH O, A1, A2, & B CELLS
- RUN AUTO CONTROL AND DAT

# ALGORITHM/GENERAL GUIDELINES FOR RESOLUTION OF ABO DISCREPANCIES

## IF PROBLEM PERSISTS:

- INCREASE INCUBATION TIME
- DECREASE TEMPERATURE OF TESTING
- RUN ANTIBODY SCREEN & PANEL, IF NECESSARY
- ENZYME TREAT CELLS, IF NEEDED
- USE OF ADSORPTION/ELUTION TECHNIQUES

---

The needs are constant.  
The gratification is instant.  
Give blood.™

• SALIVA STUDIES © D. Harming



# Where is the Discrepancy?

Anti-A Reagent	Anti-B Reagent	A1 cells Reagent	B cells Reagent	Interpret
4+	0	1+	4+	

# ALGORITHM FOR RESOLUTION

- **ALWAYS REPEAT THE TEST**
- **OBTAIN PATIENT HISTORY, AGE, DIAGNOSIS, TRANSFUSION HISTORY AND MEDICATIONS**
- **TEST PATIENT'S CELLS WITH ANTI-A,B & ANTI-A<sub>1</sub> LECTIN**
- **TEST PATIENT'S SERUM FOR ANTI-A<sub>1</sub>**
  - **Test patient serum against A1, A2, B & O cells**
  - **Run Auto Control**

# PATIENT RESULTS FROM ADDITIONAL TESTING

Patient is a 22 year-old man with no history of a transfusion.

Reagent Anti-sera				Reagent red cells				
A	B	A,B	A1 LECTIN	A1	A2	B	0	AUTO CONT
4+	0	4+	0	1+	0	4+	0	0

# INTERPRETATION

- **Patient A<sub>2</sub> with Anti-A<sub>1</sub>**
- **Transfuse Group O packed Rbc's**

**Note: Anti-A<sub>1</sub> is found in 1-8% of A<sub>2</sub> serum and in 22-35% of A<sub>2</sub>B serum. It is a cold reacting antibody and is usually clinically insignificant.**

# Main Subgroups of A: Quick Review

- **Comprises 80% of Group A Persons**
- **A<sub>1</sub>-group A cells react with both Anti-A and Anti-A<sub>1</sub>**
- **> 2 million antigenic sites per red cell**
- **Agglutinated by anti-A<sub>1</sub> Lectin**

# Main Subgroups of A

- **A<sub>2</sub> group A red cells only react with Anti-A**
- **Approx. 19 to 20% of Group A persons belong to A<sub>2</sub>**
- **Approx. 500,000 antigen sites**
- **< 1% of Group A : other weaker subgroup**
- **Red cells contain fewer antigen sites**
- **A<sub>2</sub> individuals have less transferase enzyme**
- **1 to 8% produce anti-A<sub>1</sub>**
- **22-35% produce anti-A<sub>1</sub> in A<sub>2</sub>B persons**

# MAIN SUBGROUPS OF A

Phen	Reaction of cells with reagent antiserum					Reaction of serum with reagent cells				Saliva
	A	B	A,B	H	A <sub>1</sub>	A <sub>1</sub>	A <sub>2</sub>	B	O	
A <sub>1</sub>	4+	0	4+	0	4+	0	0	4+	0	A, H*
A <sub>2</sub>	4+	0	4+	2+	0	1+/0	0	4+	0	A, H*

\* IF SE GENE IS INHERITED

# Lectins

- Proteins present in plants which bind specifically to CHO determinants and agglutinate erythrocytes by their cell surface oligosaccharide determinants
- *Dolichos biflorus*-Agglutinates A<sub>1</sub> or A<sub>1</sub>B
- *Bandeiraea simplicifolia*- B cells
- *Ulex europaeus*- H specificity

**In what order of decreasing strength  
would *Ulex europeaus*  
( *Anti-H* ) react?**

**O>A<sub>2</sub>>B>A<sub>2</sub>B>A<sub>1</sub>>A<sub>1</sub>B**

# OTHER SUBGROUPS OF A

- **A<sub>3</sub> subgroup**
  - typically demonstrates 2+mf with Anti A and Anti A,B
  - sometimes produces Anti-A<sub>1</sub>
- **A<sub>x</sub> demonstrates weak to negative reactions with Anti-A and usually 2+ reactions with Anti-A,B; usually makes Anti-A<sub>1</sub>**
- **All other weaker subgroup, A specificity can only be demonstrated by absorption/elution procedures**

# OTHER SUBGROUPS OF A

Phen	Reaction of cells with reagent antiserum					Reaction of serum with reagent cells				Saliva
	A	B	A,B	H	A <sub>1</sub>	A <sub>1</sub>	A <sub>2</sub>	B	O	
A <sub>3</sub>	2+ mf	0	2+ mf	3+	0	1+ / 0	0	4+	0	A, H*
A <sub>m</sub>	0 / +/-	0	0 / +/-	4+	0	0	0	4+	0	A, H*
A <sub>x</sub>	0 / +/-	0	1+ / 2+	4+	0	2+ / 0	0	4+	0	H
A <sub>el</sub>	0	0	0	4+	0	2+ / 0	0	4+	0	H

IF SE GENE IS INHERITED

The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening



# Subgroups of B: A Quick Review

Phen	Reaction of cells with reagent antiserum				Reaction of serum with reagent cells				Saliva
	A	B	A,B	H	A <sub>1</sub>	A <sub>2</sub>	B	O	
B	0	4+	4+	2+	4+	4+	0	0	B,H*
B <sub>3</sub>	0	1+ mf	2+ mf	3-4+	4+	4+	0	0	B,H*
B <sub>m</sub>	0	0	0/ +/-	3-4+	4+	4+	0	0	B,H*
B <sub>x</sub>	0	0/ +/-	0/ 1+	3-4+	4+	4+	0	0	H

IF SE GENE IS INHERITED

The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening



# Where is the Discrepancy?

Anti-A Reagent	Anti-B Reagent	A1 cells Reagent	B cells Reagent	Interpret
0	4+	W+	0	

# RESOLUTION

- **ALWAYS REPEAT THE TEST**
- **OBTAIN PATIENT HISTORY, AGE, DIAGNOSIS, TRANSFUSION HX & MEDICATIONS**
- **INCUBATE PATIENT'S PLASMA AND REAGENT CELLS FOR 15-30 MINUTES AT ROOM TEMPERATURE OR 15 MIN. AT 4<sup>0</sup> C**
- **ALWAYS RUN O CELLS AND AUTO CONTROL**

# PATIENT RESULTS FROM ADDITIONAL TESTING

**Eighty-six year old patient with a bleeding ulcer and no history of transfusion**

## RESULTS OF RT INCUBATION

### FORWARD

### REVERSE

Anti-A	Anti-B	A <sub>1</sub> cells	B cells	O cells	Auto contro
0	4+	3+	0	0	0

- Group B
- Elderly patient  $> 80$  years old
- Transfuse Group B Rbc's, if Antibody screen is negative

# Where is the Discrepancy?

<b>Anti-A Reagent</b>	<b>Anti-B Reagent</b>	<b>A1 cell Reagent</b>	<b>B cell Reagent</b>	<b>Interpret</b>
<b>0</b>	<b>4+</b>	<b>4+</b>	<b>2+</b>	

# Unexpected Rxns (Reverse)

- Cold Reacting Alloantibody  
(i.e. Anti-M, Anti-P<sub>1</sub> most common)
- Cold Reacting Autoantibody  
(i.e. Anti-I, Anti-H, Anti-IH)
- Passively Acquired Antibody  
(i.e. plasma exchange, mismatched Platelets)

# ADDITIONAL TESTING AND INVESTIGATION

- ALWAYS REPEAT THE TEST
- OBTAIN PATIENT HISTORY, AGE, DIAGNOSIS, TRANSFUSION HX & MEDICATIONS
- RUN O CELLS, AUTO CONTROL, AB SCREEN

FORWARD			REVERSE			
A	B	A,B	A1	B	O	AUTO CONT
0	4+	4+	4+	2+	2+	0

The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening

# ANTIBODY SCREENING RESULTS

- **POSITIVE ANTIBODY SCREEN**
- **NEGATIVE AUTO CONTROL**

## **THINK: ALLOANTIBODY**

- **Perform antibody ID**
- **Type reagent B cells for the specific antigen of the identified antibody to explain the reagent B cell reaction in the reverse grouping**

---

## **IF AUTO CONTROL IS POSITIVE**

### **THINK: AUTOANTIBODY**

- **Perform cold panel, autoabsorption if patient has **not** been transfused within the last 3 months, or alloabsorption using REST**
- **Repeat Reverse grouping using the absorbed serum or repeat reverse testing at 37<sup>0</sup>**
- **Run panel on absorbed serum to detect any underlying cold or RT reacting alloantibodies**

# Where is the Discrepancy?

Anti-A Reagent	Anti-B Reagent	A1 cell Reagent	B cell Reagent	Interpret
4+	4+	2+	2+	

**If the patient is AB Rh pos, you would need to run a saline control.**

- **ALWAYS REPEAT THE TEST**
- **OBTAIN PATIENT HISTORY, AGE, DIAGNOSIS, TRANSFUSION HX & MEDICATIONS**
- **RUN ANTIBODY SCREEN & PANEL IF NECESSARY**

# REPEAT TESTING

Reaction of red cells with reagent antiserum			Reaction of serum with reagent cells			
A	B	A,B	A <sub>1</sub>	B	O	Auto ct
4+	4+	4+	2+	2+	2+	2+

**IS THIS:**

**Group AB--Cold autoantibody (anti-I)?**

**Group AB--Cold autoantibody (anti-I) and cold alloantibody (anti-M, P<sub>1</sub>, Le<sup>a</sup>, Le<sup>b</sup>) ?**

**Group AB—Rouleaux?**

# PATIENT HISTORY

**PATIENT IS AN 82 YEAR OLD BLACK MAN WITH A DIAGNOSIS OF MULTIPLE MYELOMA WITH NO HISTORY OF A TRANSFUSION**

**CONSIDER:**

**ROULEAUX FORMATION:** “**stack of coins**” appearance of agglutination under the microscope

- High protein concentration in patient serum alters net negative charge on RBC yielding pseudoagglutination

# ADDITIONAL TESTING AND INVESTIGATION

- **Perform Saline Replacement** (assuming cell suspensions are not routinely washed ).
  - Remove serum from test tube & replace with equal number of drops of saline.

Results of saline replacement

<b>ANTI-A</b>	<b>ANTI-B</b>	<b>A1 CELLS</b>	<b>B CELLS</b>	<b>AUTO</b>
<b>4+</b>	<b>4+</b>	<b>0</b>	<b>0</b>	<b>0</b>

# INTERPRETATION

- GROUP AB
- TRANSFUSE GROUP AB

- **Other Plasma Cell Dyscrasias**
- **Wharton's jelly ( newborn Cord samples)**
- **Dextran or hydroxy ethyl starch IV infusion**  
**(Crosslinks RBCs yielding**  
**pseudoagglutination)**

# Where is the Discrepancy?

Anti-A Reagent	Anti-B Reagent	A1 cell Reagent	B cell Reagent	Interpret
0	0	0	0	

- **PATIENT HISTORY IS IMPORTANT**

# RESOLUTION

## PATIENT HISTORY

- **NEWBORN REQUIRING SURGERY FOR A HEART DEFECT**
- **GROUP O**

**NOTE: Neonatal samples**

**– Serum testing not commonly performed**

# Unexpected Rxns (Reverse)

- **Elderly**
- **Newborns**
- **Hypogammaglobulinemia (e.g. CLL, malignant lymphoma,**
- **Immunodeficiency diseases**
- **Transplant Patients**
- **ABO Subgroups**
- **Patients who received plasma transfusions or exchanges**

# Where is the Discrepancy?

Anti-A Reagent	Anti-B Reagent	A1 cells Reagent	B cells Reagent	Interpret
2+mf	0	0	4+	

The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening

# ABO DISCREPANCY: RED CELL

## PROBLEM

### CAUSE OF DISCREPANCY?

- Group A recently transfused with Group O red cells
- Group A recently transplanted with Group O BM or PBSC
- Group A<sub>3</sub> subgroup that exhibits characteristic mixed-field agglutination

Anti-A Reagent	Anti-B Reagent	A1 cells Reagent	B cells Reagent	Interpret
2+mf	0	0	4+	

# RESOLUTION

- **ALWAYS WASH PT CELLS, REPEAT THE TEST**
- **OBTAIN PATIENT HISTORY, AGE, DIAGNOSIS, TRANSFUSION HX & MEDICATIONS**
- **TEST PATIENT'S CELLS WITH LECTINS IF AVAILABLE**
- **RUN DAT AND AUTO**

# PATIENT RESULTS FROM ADDITIONAL TESTING

The patient is a 16 year-old boy with no history of a transfusion.

Reagent Anti-sera				Reagent red cells				
A	B	A,B	H LECTIN	A1	B	0	DAT	AUTO CONT
2+mf	0	2+mf	3+	0	4+	0	0	0

# INTERPRETATION

- **A<sub>3</sub> Subgroup**
- **Transfuse Group O packed Rbc's**

**This case shows the importance of transfusion history.**

# Where is the Discrepancy?

<b>Anti-A Reagent</b>	<b>Anti-B Reagent</b>	<b>A1 cell Reagent</b>	<b>B cell Reagent</b>	<b>Interpret</b>
<b>4+</b>	<b>1+</b>	<b>0</b>	<b>4+</b>	

# RESOLUTION

- ALWAYS WASH PT CELLS, REPEAT THE TEST
- OBTAIN PATIENT HISTORY, AGE, DIAGNOSIS, TRANSFUSION HX & MEDICATIONS
- RUN AUTOCONTROL
- TEST CELLS WITH HUMAN DERIVED ANTI-B REAGENT THAT HAS BEEN ACIDIFIED TO PH OF 6.0 (ACIDIFIED ANTI-B REACTS ONLY WITH TRUE “B” ANTIGEN)
- OR
- TEST PT CELLS WITH *BANDEIRA SIMPLICIFOLIA* (BS-1) IF AVAILABLE (LECTIN ONLY REACTS WITH TRUE “B” ANTIGEN)

# PATIENT RESULTS FROM ADDITIONAL TESTING

**PATIENT HISTORY: PATIENT HAS INTESTINAL BLOCKAGE AND REQUIRES SURGERY**

<b>ANTI-A</b>	<b>ANTI-B</b>	<b>A1 CELLS</b>	<b>B CELLS</b>	<b>AUTO CONT</b>	<b>BS-1 lectin</b>
<b>4+</b>	<b>0</b>	<b>0</b>	<b>4+</b>	<b>0</b>	<b>0</b>

- NOTE: Patient's anti-B will NOT react with: acquired "B" and BS-1 – lectin reacts with true "B" antigen

# INTERPRETATION

- **GROUP A with Acquired B Phenomenon**
  - Intestinal blockage→ back-up of E. coli
  - E. coli enzyme modifies “A” antigen into “B” like specificity
  - Seen in Group A individuals
- **Acriflavin dye**
  - A few patients, on rare occasions, have antibodies against acriflavin, a yellow dye used in some commercial anti-B reagent. The acriflavin-antiacriflavin complex attaches to the patient’s RBCs causing agglutination in the forward testing.

# Unexpected Rxns (Forward)

- **Out of Group Transfusion (i.e. Group O transfused to an A or B patient)**
- **Out of Group Transfusion Bone Marrow Transplantation**
- **Leukemia/Lymphoma**
- **Subgroups**
- **Hodgkins Disease**
- **Acquired B Phenomena**
- **Warm autoantibodies**

# Where is the Discrepancy?

<b>Anti-A Reagent</b>	<b>Anti-B Reagent</b>	<b>A1 cell Reagent</b>	<b>B cell Reagent</b>	<b>Interpret</b>
<b>0</b>	<b>0</b>	<b>4+</b>	<b>0</b>	

# RESOLUTION

- **ALWAYS WASH PT CELLS**
- **REPEAT THE TEST**
- **OBTAIN PATIENT HISTORY, AGE, DIAGNOSIS, TRANSFUSION HX & MEDICATIONS**

# PATIENT RESULTS FROM ADDITIONAL TESTING

PATIENT IS A 26 YEAR OLD MALE WITH ADENOCARCINOMA IN THE LIVER, STOMACH, AND INTESTINES.

## WASHED PT RBCS

ANTI-A	ANTI-B	A1 CELLS	B CELLS	INT	
0	3+	4+	0		

# INTERPRETATION

- **Group B, transfuse Group B blood**
- **Certain tumors release increased amounts of soluble substance with A &/or B substance**  
**(Adenocarcinomas of pancreas, stomach, ovary & biliary system)**
- **Soluble A or B substances neutralizes anti-A or Anti-B reagent in Forward Grouping**

# Where is the Discrepancy?

Anti-A Reagent	Anti-B Reagent	A1 cell Reagent	B cell Reagent	Interpret
0	4+mf	3+	0	

## POSSIBILITIES?

- GROUP B PREVIOUSLY TRANSFUSED WITH GROUP O
- GROUP B WITH O BMT
- GROUP B EXCHANGE TRANSFUSION WITH GROUP O
- B<sub>3</sub> SUBGROUP

- **ALWAYS WASH PT CELLS**
- **REPEAT THE TEST**
- **OBTAIN PATIENT HISTORY, AGE, DIAGNOSIS, TRANSFUSION HX & MEDICATIONS**
- **RUN AUTO CONTROL, O CELLS, DAT**

# RESOLUTION

- PATIENT HISTORY:**

- PATIENT IS A 45 YEAR OLD FEMALE ADMITTED TO SHOCK TRAUMA WITH MASSIVE INTERNAL INJURIES FROM A CAR ACCIDENT**
- SIX UNITS OF O NEG BLOOD AND 2 UNITS OF FFP FROM B+ DONOR WAS GIVEN TO THE PATIENT**

## FORWARD

## REVERSE

A	B	A,B	H LECTIN	A1	B	0	DAT	AUTO CONT
0	4+mf	4+mf	3+	4+	0	0	0	0

# INTERPRETATION

- **TRANSFUSION OF GROUP O RBCs TO B PATIENT**
- **THE MOST COMMON CAUSE OF MIXED FIELD AGGLUTINATION IS TRANSFUSION OF O CELLS TO AN A OR B PATIENT**

# OTHER CAUSES OF MIXED FIELD AGGLUTINATION

- **OUT OF GROUP BONE MARROW TRANSPLANTATION**
- **FETAL-MATERNAL BLEED**

# UNEXPECTED RESULTS IN THE FORWARD AND REVERSE GROUPING

- **CHIMERISM**



The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening



# Twins

	Anti-A Reagent	Anti-B Reagent	Anti- A,B	A1 cell	B cell
Twin1	0	2+mf	2+mf	4+	0
Twin 2	0	+wk	+wk	4+	0

- Twin 1      70% B   30% O
- Twin 2      30% B   70% O

# NETWORKING AND COOPERATION GIVES SOLUTIONS TO PROBLEMS



The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening


# THE END



The need is constant.  
The gratification is instant.  
Give blood.™

© D. Harmening

# INHERITANCE OF Cis-AB



**Figure 5-19.** Example of cis-AB inheritance to unequal crossing-over. □ = male; ○ = female. (From Harmening-Pittiglio,<sup>33</sup> p 7 with permission.)